

Table 1. The immunophenotyping panel used for flow cytometry results presented in Figure 1.

Marker	Fluorophore	Clone	Vendor
CD45	PerCP	30-F11	BD
B220	APC	RA3-6B2	eBioscience
CD19	APCeFl780	1D3	eBioscience
IgM	PECF594	R6-60.2	BD
IgD	eFl450	11-26	eBioscience
CD8	BV510	53-6.7	BioLegend
CD4	PECy5.5	GK1.5	SouthernBiotech
CD44	BV785	IM7	BioLegend
CD62L	BV570	MEL-14	BioLegend
CD25	BB515	PC61	BD
F4/80	PECy7	BM8	eBioscience
CD11b	PerCPCy5.5	M1/70	BD
CD11c	AF647	N418	BioLegend
CD138	PE	281-2	BD
Ly6c	BV711	HK1.4	BioLegend
NK1.1	BV480	PK136	BD
CD5	BV605	53-7.3	BD
CD21	FITC	4E3	eBioscience
CD23	BUV737	B3B4	BD
Zombie NIR Fixable Viability Dye			BioLegend